# **CAXS MEDICAL**

PRESENTS

# **3D MODELIZATION**



BIOMOD 3S BIOMOD MI AND BIOMOD FB



## **BIOMOD 3S**

**THE BIOMOD 3S SYSTEM IS AN INNOVATIVE OPTICAL IMAGING TECHNOLOGY** combined with an X-ray system that provides a 3D modelization of the spine with patients in a natural standing position.

BIOMOD 3S is used for the evaluation and diagnosis of deforming pathologies of the spine and balance disorders in children and adults.

This solution is compatible with rotational single-source radiology systems.



#### **EXAMINATION**

The examination of a patient is performed in a functional standing position. Two optical acquisitions of the back (only takes a few seconds), are coupled with frontal and sagittal radiographs (stitching, frontal & profile) to obtain the 3D model. Four radio-opaque markers are positioned on the back (C7, L5 and EIPS).The vertebral column is drawn with the red marker.

#### ACQUISITION

The optical acquisition of the BIOMOD column is instantaneous and performed during the two radiographic acquisitions (frontal/ sagittal). The optical and radiographic acquisitions of the frontal/profile images are merged to create the 3D model of the column. There is no change to the classic clinical routine for the practitioner or the patient.

## RECONSTRUCTION

The station merges optical acquisitions with radiological acquisitions in a few steps, to obtain a 3D model of the spine and back. The steps include:

1) Identification of markers and the position of the spinal line on the optical acquisition

2) Tracing the position of the spine on radiographic acquisition

3) Placement and adjustment of the vertebrae

The reconstruction time is about 10 minutes.

#### **ADDED VALUE**

• 2D frontal and profile parameters (provides additional information on twists and vertebral rotations compared to conventional 2D parameters measured on radiographs).







 Biomod 3S technology can be integrated into radiology rooms equipped with a system that produces full-spine digital images.



## **BIOMOD MI**

#### BIOMOD MI IS A MEDICAL SYSTEM THAT OFFERS A 3D MODEL OF LOWER LIMBS IN A FUNCTIONAL STANDING POSITION.

The non-invasive optical BIOMOD technology, couples two standard frontal and sagittal radiographs with optical acquisitions of the lower limbs to obtain a 3D model.

BIOMOD MI is a module of the BIOMOD platform that also includes the FB and the 3S module.



#### EXAMINATION

The Biomod MI provides access to a 3D visualization of bone structure and offers a complete 3D evaluation of the lower limbs in a functional upright position. It provides a better understanding of bone misalignment and rotations/torsions. With the Biomod MI, it is possible to obtain better accuracy\* of the automatically calculated clinical parameters, and save time\* throughout the patient's journey, from diagnosis to postoperative follow-up.

#### **ACQUISITION**

Optical images X-ray images



#### RECONSTRUCTION

The software merges optical acquisitions with radiological acquisitions in a few steps to obtain a 3D model of the lower limbs according to the following characteristics: radio-opaque marker tracking, anatomical landmark identification and clinical exploitation.

#### **CLINICAL PARAMETERS**

	<ul> <li>Femoral head diameter</li> <li>Femoral axis length</li> <li>Length tibia axis</li> <li>Total length</li> </ul>
TORSIONS / ROTATIONS	<ul><li>Femoral torsion</li><li>Femoro-tibial rotation</li></ul>
	<ul> <li>Varus/Valgus</li> <li>Flessum/Recurvatum</li> <li>Angle HKS</li> <li>Mechanical femoral angle</li> <li>Anatomical femoral angle</li> <li>Mechanical tibial angle</li> </ul>



# **BIOMOD FB**

#### BIOMOD FB IS A DEVICE FOR ANALYZING POSTURE AND BALANCE.

This solution provides almost instant access to the clinical parameters used to monitor dorsopathies, characterization of lower limb deformities and posture disorders.

#### **FEATURES**

- 3D characterization of the entire body
- Objective measurements and calculations of clinical parameters
- Morphometric report
- Comparison and history of changes in clinical parameters
- Educational and interactive, documents for the patient file



#### **CLINICAL PARAMETERS**

- Follow-up of scoliosis / kyphosis, sagittal balance disorders and degenerative spine
- Monitoring and study of standing posture
- Evaluation and monitoring of rehabilitation programs and orthopaedic treatments
- Pre- and post-processing comparison tool
- Visualization of correction effects on the patient's morphology





#### SETTINGS

- Objective quantitative measurements of static and posture under load
- Whole body mapping (level lines and concavities)
- Deviations and imbalances (pelvis, shoulders)
- Back and lumbar arrows
- Deformations of the lower limbs



## CONTRIBUTION

- Contribution to the reduction of patients' exposure to medical irradiation (up to -50% use of ionizing radiation)
- Assistance in diagnosis and decision-making (quantified results and recognized data)
- Monitoring the effectiveness of treatment (visualization of the evolution, patient adherence)
- Monitoring of rehabilitation programs (patient file with quantified and verifiable progress)
- Improved communication with the patient (visualization of the patient's back, treatment results for complete adherence to the protocol)
- Modern assessment tool (innovative technological contribution to clinical practice. Minimal physical contact with the patient)

#### 3D MODELIZATION : MADE IN FRANCE INNOVATION



HEADQUATERS > 393 rue Charles Lindbergh – 34130 Mauguio – FRANCE FACTORY > Parc Scientifique G. Besse – 175 Allée Von Neumann – 30035 Nîmes Cedex 01 – FRANCE PHONE > +33 4 66 29 09 07 – FAX > +33 4 67 50 49 09 – www.dms-imaging.com

These products are Class IIa medical devices manufactured by AXS MEDICAL and bears the CE marking. They are certified in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, exclusively for the indication of balance and posture analysis. Other non-medical uses ascribed to these devices are not within the scope of CE certification, and user should be aware product performance and/or safety has not been evaluated by SGS for those purposes. Read the operating instructions carefully. Photos are non-contractual. BrochureA 3D MODELIZATION 02 05/2020